

FloodLAMP

Biotechnologies, PBC


Extending the Reach of LAMP:
the FloodLAMP platform for globally accessible disease screening

New England Biolabs Seminar – March 3, 2022

Randy True, Founder and CEO | randy@floodlamp.bio



David Deutsch's Deep Optimism:

a much more general truth, namely that, for people, *problems are inevitable*. So let us carve *that* in stone:



It is inevitable that we face problems, but no particular problem is inevitable. We survive, and thrive, by solving each problem as it comes up. And, since the human ability to transform nature is limited only by the laws of physics, none of the endless stream of problems will ever constitute an impassable barrier. So a complementary and equally important truth about people and the physical world is that

problems are soluble. By 'soluble' I mean that the right knowledge would solve them. It is not, of course, that we can possess knowledge just by wishing for it; but it is in principle accessible to us. So let us carve that in stone too:



The success of the Tech giants — Amazon, Apple, Microsoft, Alphabet, Facebook, Google, etc. — confined to their homes, they turned to the internet for pretty much everything:

The Coronavirus
Pandemic >

LIVE Covid-19 Updates

All Anyone Wants for Christmas Is a Covid Test

Bloomberg
Businessweek

Want to go to a restaurant?

We need better testing

A sports event? A movie?

We need better testing

Travel somewhere? Anywhere?

We need better testing

Go to a wedding? A funeral? A party?

We need better testing

Send your kids back to an actual school?

We need better testing

A normal-ish life before there's a vaccine?

We need better testing

Covid has created a capitalist nightmare

Big Tech and Big Pharma are more powerful
than ever

The New York Times


Over the past two years, as the pandemic claimed the lives of millions of

people and upturned the lives of everyone else, a silent revolution was taking place. Western capitalism suddenly found itself usurped; replaced by an even more oligarchic and authoritarian capitalist mode of power — what we might call neo-feudalism. Even as restrictions are lifted in the UK and other countries, the consequences of these changes will be felt for a long time to come.

On the one hand, hundreds of thousands of small and medium businesses were forced to shut down, triggering one of the worst jobs crises since the Great Depression. In 2020 alone, workers around the world cumulatively lost \$3.7 trillion in earnings — an 8.3% decline.

**Demand for Covid-19 testing is falling, but experts
caution it's as important as ever**

“FloodLAMP's mission is to improve global health and resiliency through universal access to rapid molecular testing.”



1) Real World Programs

Operating Programs in the Real World



EMS Leadership Conference Screening

→ 1,000 attendees



Municipal & EMS Deployments

- Coral Springs, FL approx. 132,000
- Davie, FL approx. 104,000
- Bend, OR approx. 94,000




Preschool Screening with Family Pooling


- Approx. 40 families
- 65-105 people per bi-weekly screening




Town of Davie, FL



Coral Springs, FL



Coral Springs, FL




↙ This novice (but dedicated) tester is able to screen 1000 critical city employees.

Davie, FL




Test Diagrams

FloodLAMP QuickColor™ COVID-19 Test





Bend, OR



Nova, Southeastern University Lab

Lab setup (Jan 2022)



FloodLAMP accepted into Levan Innovation Center



Early Childcare COVID-19 Family Screening Pilot: Carillon Preschool

Protecting young children during the Omicron surge

Prepared by:

Randy True – CEO, FloodLAMP Biotechnologies, PBC
randy@floodlamp.bio

Theresa Ling – User Experience and Design, FloodLAMP Biotechnologies, PBC
theresa@floodlamp.bio



"I really appreciate this process. D tested positive on antigen test at home today and is now showing symptoms, so I truly appreciate the surveillance program picking this up! Thank you again!!"

- Parent of a student who was referred to follow up testing while negative on antigen and pre-symptomatic

"Thank you so much for helping to get our kids tested and doing what you can to keep them safely in school! We are all so lucky to have you as part of our community!"

- Parent of a student and MD, wrote to school administrators to move to mandatory testing with FloodLAMP

"I thought [FloodLAMP] was perfect. Easy and convenient."

- Parent of a student and restaurant owner who later wanted to extend FloodLAMP testing to her restaurant employees


"You are a godsend to our community. I cannot begin to thank you for all of your support and love in keeping our community safe in this difficult time."

- Preschool teacher and administrator

[Link for Preschool Pilot Description](https://docs.google.com/document/d/1TOnklq-65XUX-v-li018rteK9jeL8NgqbNrtbndgD_o)
https://docs.google.com/document/d/1TOnklq-65XUX-v-li018rteK9jeL8NgqbNrtbndgD_o

FLOODLAMP ARCHIVE FILE PATH:
various/fl-whitepapers/FloodLAMP Whitepaper -
California Preschool Family Pooled Screening
Pilot (June 2022).md

Preschool, CA




Pre-School Screening – Case Report

- Several infections detected by FloodLAMP 1-2 days before symptoms or antigen tests
- In 3 different cases, family members of the preschool student were detected, showing the advantage of family pooling in protecting spread in a congregate setting such as schools and workplaces.

	1/2/22	1/3/22	1/4/22	1/6/22	1/7/22	1/8/22	1/9/22	1/10/22	1/13/22	1/14/22	1/17/22
	0	1	2	4	5	6	7	8	11	12	15
Family Pool	POS FloodLAMP										POS FloodLAMP
Parent 1	Asymptomatic										
	POS FloodLAMP (Family Pool)	POS FloodLAMP (Indiv)					POS FloodLAMP (Indiv)	> POS FL Result			POS FloodLAMP (Family Pool)
	NEG Antigen									NEG Antigen	
			POS PCR Lab >	> > > >	> POS PCR Result						
Parent 2	Asymptomatic		Symptoms								
	POS FloodLAMP (Family Pool)	POS FloodLAMP					POS FloodLAMP (Indiv)	> POS FL Result			POS FloodLAMP (Family Pool)
	NEG Antigen									NEG Antigen	
		POS PCR Lab >	> POS PCR Result								
Child 1	Asymptomatic		Symptoms	Asymptomatic							
	POS FloodLAMP (Family Pool)	NEG FloodLAMP					POS FloodLAMP (Indiv)	> POS FL Result			POS FloodLAMP (Family Pool)
			POS Antigen							NEG Antigen	
					POS PCR Lab >	> POS PCR Result					
Child 2	Asymptomatic		Symptoms	Asymptomatic							
	POS FloodLAMP (Family Pool)	NEG FloodLAMP					POS FloodLAMP (Indiv)	> POS FL Result			POS FloodLAMP (Family Pool)
			POS Antigen						POS Antigen (slight)	NEG Antigen	
				POS PCR Lab >	> POS PCR Result						

Pre-School Screening



Mon 1-10



Thurs 1-13

Mix of self collected individuals and family pools

New Year's Eve Case Report

- FloodLAMP caught asymptomatic infection of advisor that 3 expensive OTC/POC rapid molecular tests missed (Detect, Lucira, and Mesa Accula)
- BinaX now wasn't positive until 2 days later.
- All FloodLAMP positive samples confirmed by PCR


FloodLAMP New Years Eve Case Report caught asymptomatic infection that all antigen and OTC/POC molecular tests missed						
FloodLAMP EasyPCR on saliva				POS - SalivaDirect PCR Ct 27.7 8 PM	POS - SalivaDirect PCR Ct 29.0 11 AM	
Accula Rapid Molecular PCR (Worksite)					Accula NEG - nasal swab 4 PM	
Lucira Rapid Molecular LAMP				Lucira NEG - nasal swab 5 PM		
Detect Rapid Molecular LAMP	NEG - nasal swab 10 AM		Detect NEG - nasal swab 2 PM			
CaseStart Antigen						
FlowFlex Antigen		NEG - nasal swab 8 AM	NEG - nasal swab 2 PM			POS - nasal swab 10 AM
BinaXNow Antigen				NEG - nasal swab 5 PM	NEG - nasal and throat swabs 10 AM	POS - nasal swab 10 AM
FloodLAMP QuickColor LAMP Nasal Swab		NEG - nasal swab 1 PM		NEG - nasal swab PCR Ct Und/39.4 5PM	POS - nasal swab PCR Ct 32.2 10 AM	
FloodLAMP QuickColor LAMP Throat Swab			POS - nasal swab PCR Ct 33.9 9 AM	POS - throat swab PCR Ct 30.2 5 PM	POS - throat swab PCR Ct 35.5 10 AM	
	12-27-21	12-29-21	12-31-21	12-31-22	1-1-22	1-2-22

Surveillance Testing Stats

FloodLAMP Surveillance Testing Stats			Thru Mar 2, 2022
<u>Org</u>	<u>Pools</u>	<u>People</u>	<u>Positives</u>
FL FF	1,107	1,779	44
FL FTFC	61	195	0
Kent Summer Camp	190	696	0
Coral Springs City	7,400	21,676	348
Davie Fire Rescue	2,235	4,226	55
Pink Shoes Production	671	671	1
Bend Fire Rescue	617	617	215
TOTAL	12,281	29,860	663
Avg Pool Size	2.4		

2) Platform – Wrap Around Components


2) Platform – Wraparound program components



Reagent Kits



Collection Kits



Collection Instructions



**FloodLAMP
Biotechnologies**
A Public Benefit Corporation

FloodLAMP Pooled Swab
Collection Kit

For use with the FloodLAMP
Mobile App



Open your phone camera
and scan to see these
instructions online

FloodLAMP.bio/collectionkits

1

Begin collection
Sponsors: Tap the "Collect" button from the app home screen. Follow the instructions in the "Collect your samples" screen, which are listed here below.

Participants: Follow the swabbing instructions listed below in coordination with your group's sponsor.

2

Sanitize
Wash and dry your hands before collecting your sample, or use hand sanitizer.

3

Open the kit
Sponsors: Open your kit and distribute swabs to each participant who will be adding samples to the tube.

4

Open swab packs
Each participant should open their swab pack. Be careful not to touch the absorbent end with your hands.

5

Swab the first nostril
Each participant should insert the swab tip into their first nostril until the absorbent end is no longer visible. Swab the sides of the nostril in a circular motion 4 times. Make sure to twist the swab as you go.

6

Swab the other nostril
Each participant should repeat the previous steps with the same swab in the second nostril.

7

Put swabs in tube
Each participant should snap off the end of the swab at the break point. Be careful not to touch the absorbent end with your hands.

Drop the shortened swab into the barcoded tube.

8

Cap tube
Sponsors: Once each participant has added their swab, securely tighten the cap onto the tube. Tap "Next" at the bottom of the "Collect your samples" screen.

9

Scan tube
Sponsors: On the "Scan" screen, capture the code on the side of the collection tube. Tap "Next" to continue.

10

Add participants
Sponsors: On the "Add Participants" screen, add the names of the pool participants (including the Sponsor if they swabbed) by searching for their names or FloodLAMP registered email addresses. You can add up to 4 swabs in each tube. Participating minors will be registered under their guardian's account.

If you have trouble locating a participant in the system, add their name in the "Collection Notes" field and make sure they sign a consent form.

11

Put the tube in the biohazard bag
Securely close the zippered bag, then thoroughly wash or sanitize your hands.

12

Return the biohazard bag
By collection point: If your kit does not include a return mailer, bring your biohazard bag to the collection point designated by your program administrator.

By mail: If your kit includes a return mailer, seal the biohazard bag in the return mailer. This packaging is specially designed for carrying samples suspected of carrying infectious substances. Do NOT replace it.
Note: Your return shipping has already been paid.

13

Register your return
Sponsors: Tap the "Return" button from the app home screen.


Equipment and Supplies




"Lab"





Portable System



Portable System – in action




Prepared Cart




App

For program participants




App

For the lab




Admin portal

Tools for program administrators



Test Diagrams

FloodLAMP QuickColor™ COVID-19 Test



Training Program – Interactive Video

 FloodLAMP ENGLISH (EN) ▾

Quick Color LAMP Training

[Dashboard](#) / [My courses](#) / [QCLT](#)

Welcome

 [Welcome!](#)
To do: View

 [General Information](#)
To do: View

Interactive Video Training

 [Instructions and Expectations for Interactive Video Training](#)
Done: View

 [If you have trouble with EdPuzzle](#)
To do: View

 FloodLAMP ENGLISH (EN) ▾

5) Reaction Mix Prep

Done: View

 edpuzzle



00:00 15:36

◀ 4) Inactivation Reaction

Jump to... ▾

Training Program – Assessment for Certification

FloodLAMP
BIOTECHNOLOGIES



Video List for Phase1


- 1) Overview
- 2) Safety and Contamination
- 3) Prep 1X ISS
- 4) Prep Rxn Mix
- 5) Inactivation
- 6) Amplification (assay run)
- 7) App and Resulting
- 8) Dispenser

Moodle Course

Enrollment, completion, activity, etc data collection

EdPuzzle
Interactive
screencast w/
questions

Review Amp Run
Sheet



Quality Management System

- Implementing ISO 13485 QMS
- QR Form Driven Logs

	Document Number:	SOP-102-A_1
	Effective Date:	01/22/2022
	Version:	1
	Reaction Mix Prep - Table Form	

CLAMP MM - NEB M1804 LOT: EXP:	
<u>New In Use CLAMP MM</u> Name: Location: Date and Time:	Retire Link Name: Date and Time:

Reaction Mix Prep

- Thaw PGS and CLAMP MM
- Vortex PGS and CLAMP MM for 3s, centrifuge/shake down
- Add CLAMP MM to PGS tube (or new 1.5mL tube for custom amount)
- Vortex for 10s and centrifuge/shake down
- Add 23µL of reaction mix per well (at bottom, no blowout, no tip touch)
 - * optional use reservoir and multichannel
- Check volumes (look at pink liquid) by lifting plate or strip tubes out of rack
- Label with Reaction Mix ID

					48 8	550 92	655 110
Name	Date + Time	Reaction Mix ID (RM_MMDD-1)	CLAMP MM TubelD (_MMDD-X)	PGS LOT# + EXP	Num Rxns	PGS µL	CLAMP MM µL

Version 1

SOP-102-A Reaction Mix Prep - Table Form

Effective 22/01/2022

Confidential Do Not Distribute without Permission

Page 1 of 2

	Document Number:	SOP-103-A_1
	Effective Date:	01/22/2022
	Version:	1
	Amplification - Run Form	

PREP

- ☐ Heaters on (Heat indicator lit!)
- ☐ 1X Inactivation Saline Soln ready
- ☐ Reaction mix strip8/plates ready
- ☐ Safety Procedures:
 - lab coat
 - gloves
 - face mask
 - face shield or goggles
- ☐ Alcohol wipe sample tubes

INACTIVATION

- ☐ Add 1mL of 1X Inactiv Saline Soln to each sample tube
- ☐ Vortex 10sec (tubes) or 30sec (rack)
- ☐ Water Bath (set 99°C)
 - Heat 8 min then Cool 10 min
- OR**
- ☐ Dry Heat Block (set 95°C)
 - Heat 5 min then Cool 5 min

AMPLIFICATION REACTION

- ☐ Thaw Reaction Mix (if frozen)
- ☐ Setup tips to align with reactions
- ☐ Add 2µL from each sample tube (pipet up& down 5X, blowout in liquid, tip touch)
- ☐ Heat 65°C for 25 min, set timer
- ☐ Intake sample tubes in App
- ☐ Put sample tubes in fridge
- ☐ Remove Reactions and Let cool for 1 min before photo
- ☐ Take photo in lightbox
- ☐ Crop photo and apply vivid filter
- ☐ Update results in App
- ☐ Log Run with Form link

	Amplification Log Form Link Location:
Name:	
Date:	
Time:	

Sample Batch ID (from App):
Inactiv Heater:
1X Inactiv Saline Soln ID:
Dispense: Pipette or Manual Disp or Electric Disp

Amp Heater:	Strip8 Tubes or Plates
Reaction Mix: Frozen or Fresh	
Reaction Mix ID:	
Num Reactions (including controls):	
Pos Ctrl ID:	
Neg Ctrl ID (today's 1XISS):	
Initial Inconclusives?:	
Notes:	

Version 1


SOP-103-A Amplification - Run Form

Effective 22/01/2022

Confidential Do Not Distribute without Permission


Page 1 of 2

New Site Validation Run



3) Open Source “Generic” Tests

At Sweet Spot of Performance



Filling a Gap in the Testing Landscape



Current testing paradigms cannot stem a pandemic caused by an asymptotically transmitted, aerosolized pathogen.




- Scaling a single test up to very high levels offers **opportunities for innovation** in sampling, test chemistry, and program configuration.
- A "**Point of Need**" modality of testing is needed that combines the scalability and low cost of liquid reagent processing with the ease, flexibility and decentralized nature of POC/At-Home.
- **Public Health processing sites** can fill the gap between CLIA labs and DIY at home tests.
- Want this new capability to be **additive** and not cannibalize the clinical diagnostic infrastructure.

LAMP EUAs for SARS-CoV2

EUA Date	Test	Diagnostic	Target	IC	Detection	Extraction	Amplification	LOD spike	LOD GE/ 3ml swab
Nov-17	Lucira (At Home POC)	Lucira COVID-19 All-In-One Test Kit	N, N	External + IC	Colorimetric	~Direct	Lucira device	Inactivated virus	2,700
Oct-05	Seasun (2 duplex)	AQ-TOP COVID-19 Rapid Detection Kit PLUS	orf1ab, N	Human RNaseP	PNA Probe	Manual (Qiagen_60704 or Seasun_SS-1300) or Automated (Panagene_PNAK-1001 on PanaMax48)	BioRad_CFX96 or ABI7500	NCCP_43326 genomic RNA	3000
Sep-01	Detectachem	MobileDetect Bio BCC19 (MD-Bio BCC19) Test Kit	N, E	only external controls	Colorimetric	Direct (1µl of transport media into LAMP)	heat block or qPCR (MD-Bio heater or BioRad_T100 or AB_Veriti or ...)	Twist_MT007544.1 synthetic gRNA	225,000
Aug-31	Mammoth	SARS-CoV-2 DETECTR Reagent Kit	N	Human RNaseP	CRISPR Probe	Automated (Qiagen_955134 on Qiagen_EZ1AdvancedBenchtop)	ABI7500	SeraCare_AccuPlex_0505-0168 IVT encapsulated RNA	60,000
Aug-13	Pro-Lab	Pro-AmpRT SARS-CoV-2 Test	RdRP	only external controls	Probe	Direct (swab into 0.1ml) or kit (swab into 1ml then Pro-Lab_PLM-2000)	Optigene_GenieHT	BEI_NR-52287 inactivated virus	125*
Jul-09	UCSF/Mammoth	SARS-CoV-2 RNA DETECTR Assay	N	Human RNaseP	CRISPR Probe	automated (Qiagen_955134 on Qiagen_EZ1)	ABI7500	SeraCare_AccuPlex_0505-0126 (lot 10480311) IVT RNA	60,000
May-21	Seasun (1 duplex)	AQ-TOP COVID-19 Rapid Detection Kit	orf1ab	Human RNaseP	PNA Probe	manual (Qiagen_60704)	BioRad_CFX96 or ABI7500	NCCP_43326 gRNA	21,000
May-18/ Aug31	Color	Color Genomics SARS-CoV-2 RT-LAMP Diagnostic Assay	N, E, nsp3	Human RNaseP	Colorimetric	automated (PerkinElmer_CMG-1033 on PerkinElmer_Chemagic360)	plate reader (Biotek_NE02), Hamilton_Star	ATCC_VR-1986D (gRNA)	2,250
May-06	Sherlock BioSci	Sherlock CRISPR SARS-CoV-2 Kit	orf1ab, N	Human RNaseP	CRISPR Probe	Manual (ThermoFisher_12280050)	heat block or qPCR, Plate reader (BioTek_NE02)	gRNA	20,250
Apr-10	Atila BioSystems	LAMP COVID-19 Detection Kit	orf1ab, N	human GAPDH	Probe	Direct (swab into 350µl, 15min RT then 3µl to LAMP)	BioRad_CFX96 or ABI7500 or Roche_LightCycler480II or Atila_PG9600	SeraCare_AccuPlex_0505-0129 pseudovirus (recombinant alphavirus)	3,500*

Just in 2020!
from Matt McFarlane
(twitter @mattmcfar)

PCR and LAMP Sister Assays



Streamlined Sample Prep

- Standard Rabe-Cepko TCEP/EDTA inactivation solution
- Added directly to dry swabs
- Same sample for both tests


- Same sample for both tests
- 2 μ L sample volume



visual readout

QuickColor™ LAMP Test

- M1804 MM
- Standard by-the-book 25 μ L reaction
- 3 primer sets mixed (AS1e, N2, E1)
- 25 minutes at 65°C



Standard RT-qPCR instrument

EasyPCR™ Test


- Luna E3006 kit
- Standard by-the-book 20 μ L reaction
- CDC primers in SalivaDirect config (N1, RNaseP)

Clinical Evaluation Data

Clinical evaluation performed by the Stanford CLIA Lab, with excellent results and praise on the **"really straightforward"** protocol.

EasyPCR™ Test

- 3 copies/μl LoD
- 98% sensitivity (PPA 39/40)
- 100% Specificity (40/40)
- No false positives



Gamma inactivated cell lysate from BEI spiked into raw clinical negative sample

QuickColor™ LAMP Test


- 12 copies/μl LoD
- 90% Sensitivity (PPA 36/40)
- Missed positives only high Ct (>36 with direct PCR)
- 100% Specificity (40/40)
- No false positives

FloodLAMP SwabDirect PCR Result	Comparator Positive	Comparator Negative	Total
Positive	39	0	39
Negative	1	40	41
Invalid	0	0	0
Total	40	40	80
Positive Agreement	97.5% (39/40) 95% CI: 86.8% to 99.9%		
Negative Agreement	100% (40/40) 95% CI: 91.2% to 100%		

FloodLAMP QuickColor Test Result	Comparator Positive	Comparator Negative	Total
Positive	36	0	36
Negative	4	40	44
Total	40	40	80
Positive Agreement	90.0% (36/40) 95% CI: 76.3% to 97.2%		
Negative Agreement	100% (40/40) 95% CI: 91.2% to 100%		


Source of Specimens: Stanford COVID-19 Clinical Testing Program
Specimen Type: Anterior Nares Swab in PBS, previously tested and frozen
Comparator Test: Hologic Panther Fusion SARS-CoV-2 Assay and Hologic Panther Aptima SARS-CoV-2 Assay

Clinical Evaluation by Stanford CLIA Lab



Clinical Evaluation by Stanford CLIA Lab

Direct PCR Ct vs Colorimetric LAMP Result



40 Clinical Positive Samples:

29 Positive by FloodLAMP QuickColor™ Test


7 Inconclusive (all high Ct, low viral load)

4 Negative (all > 36 Ct, very low viral load)

40 Clinical Negative Samples:


40 Negative by FloodLAMP QuickColor™ Test

Clinical Evaluation by Stanford CLIA Lab



Clinical Evaluation by Stanford CLIA Lab

FloodLAMP Stanford Clinical Evaluation - All 3 Tests: Fluorimetric LAMP, Colorimetric LAMP, and PCR



One sample not detected by
FloodLAMP™ EasyPCR
(PPA 39/40 97.5%)

6 Positives Not Detected by
FloodLAMP QuickFluor™ test
(PPA 34/40 85%)

● Positives not detected by
FloodLAMP QuickColor™ test
(PPA 36/40 90%)
* includes initial inconclusives resolved by
PCR

All 3 Tests:
FloodLAMP EasyPCR™
QuickColor™ LAMP
QuickFluor™ LAMP
had no false positives
(NPA 40/40 100%)

Saliva Direct "most unique EUA" – FDA

- Is a protocol not a product
- Open Access - any CLIA can use
- Lots of clone assays being used

- Yale can "designate" CLIA labs

Yale SCHOOL OF MEDICINE

About YSM ▾ Education ▾ Research ▾ Patient Care ▾ Diversity ▾ News ▾ A-Z Index ▾

SalivaDirect: What You Need to Know About the New COVID-19 Test

September 03, 2020 by Michael Greenwood

COVIDTracker

About our Test Development

SalivaDirect™ has moved!

For information regarding SalivaDirect™, including lab designation and up-to-date validation results, please visit our new page at <https://publichealth.yale.edu/salivadirect/>.

COVIDTrackerCT will remain home to our surveillance, open-source protocols, and research (including test development).

Saliva

oteinase K

Combine & mix

Heat inactivate 95°C

Dualplex RT-qPCR

Created with BioRender.com

SalivaDirect: Lab authorization request form

SalivaDirect received Emergency Use Authorization (EUA) from the Food and Drug Administration (FDA) on August 15th, 2020.

The SalivaDirect FDA EUA is for a Laboratory Developed Test, but we have the right to designate other labs its use.

Only high complexity CLIA-certified labs in the United States can become authorized to run SalivaDirect. Please fill out this form if you would like to receive more information on how to become a designated lab.

Do you have access to at least 5x SARS-CoV-2 positive and 5x negative saliva samples for CLIA verification? *

Would you be open to collaboration on bridging studies for your automation systems?

Open EUA's – 2 Full Submissions + Pre EUA

FloodLAMP QuickColor™ COVID-19 Test

Instructions for Use v1.2

IVD

COVID-19 Emergency Use Authorization Only
For *in vitro* diagnostic (IVD) Use

www.floodlamp.bio
FloodLAMP Biotechnologies, PBC | 930 Brittan Ave. San Carlos, CA 94070 USA

FloodLAMP EasyPCR™ COVID-19 Test

Instructions for Use v1.1

IVD

COVID-19 Emergency Use Authorization Only
For *in vitro* diagnostic (IVD) Use

www.floodlamp.bio
FloodLAMP Biotechnologies, PBC | 930 Brittan Ave. San Carlos, CA 94070 USA



**FloodLAMP
Biotechnologies**

A Public Benefit Corporation

**FloodLAMP Pooled Swab
Collection Kit DTC**

For use with the FloodLAMP
Mobile App

Open EUA – Disclosure of Components

FloodLAMP QuickColor™ COVID-19 Test

IVD FDA Emergency Use Authorization - Instructions for Use (v1.2)

Table 1: Validated reagents used with the Test

Item	Concentration	Chemical Composition	Vendor	Catalog Number
TCEP	.5 M	tris(2-carboxyethyl)phosphine hydrochloride	Sigma-Aldrich	646547-10X1ML
EDTA	.5 M	Ethylenediaminetetraacetic acid	Thermo Fisher	15575020
NaOH	10 N	Sodium Hydroxide	Sigma-Aldrich	SX0607N-6
Nuclease-free Water		Ultrapure Water, nuclease-free	Thermo Fisher	10977015
NaCl	5 M	Sodium Chloride	Thermo Fisher	24740011
Guanidine HCl	6 M	Guanidine Hydrochloride	Sigma-Aldrich	SRE0066
Colorimetric LAMP MM*		Colorimetric LAMP Master Mix	New England Biolabs	M1804

* Item may not be substituted for equivalent. Only the specified vendor and catalog number may be utilized.

The FloodLAMP QuickColor™ COVID-19 Test uses 18 LAMP primers targeted for 3 different SARS-CoV-2 genes, with 6 primers for each target. All 18 primers are mixed together and are input into a single amplification reaction. Primer names and sequences are listed in Table 2. Primers may be purchased pre-blended from the vendor LGC Biosearch Technologies with the product names LAMP_S2-As1e, LAMP_S2-N2, LAMP_S2-E1. Alternatively, primers may be purchased as individual custom oligos. Appropriate validation of primer mixes from custom oligos is required. See Primer Preparation below for more information.

Table 2: Primer names and sequences

Primer Name	Sequence (5'-3')
ORF1ab gene (AS1e)	
Orflab_FIP	TCAGCACACAAAGCCAAAATTTATTTTCTGTGCAAAGGAAATTAAGGAG
Orflab_BIP	TATTGGTGGAGCTAAACTTAAAGCCTTTCTGTACAATCCCTTTGAGTG
Orflab_F3	CGGTGGACAAATTGTAC
Orflab_B3	CTTCTCTGGATTAAACACACTT
Orflab_LF	TTACAAGCTTAAAGAATGTCTGAACACT
Orflab_LB	TTGAATTTAGGTGAACATTGTGCAG
N Gene (N2)	
N2_FIP	TTCCGAAGAACGCTGAAGCGGAAGTATTACAAACATTGGCC
N2_BIP	CGCATTGGCATGGAAGTCACAAATTTGATGGCACCTGTGTA
N2_F3	ACCAGGAACATACAGACAA
N2_B3	GACTTGATCTTTGAAATTTGGATCT
N2_LF	GGGGGCAAAATTGTGCAATTTG
N2_LB	CTTCGGGAACGTGGTTGACC

Table 7: Primer-Guanidine Solution

Component	Volume (1 reaction)	Volume (1 reaction x 100) 1 x 96-plate w/ 4% overage
10X LAMP Primer Mix	2.5 µL	250 µL
Guanidine HCl (400 mM)	2.5 µL	
Guanidine HCl (6 M)		16.7 µL
Nuclease-free Water	5.5 µL	783 µL
TOTAL VOLUME	10.5 µL	1050 µL

Table 8: Colorimetric LAMP Amplification Reaction

Component	Volume (1 reaction)	Volume (100 reactions)
Primer-Guanidine Solution	10.5 µL	1050 µL
Colorimetric LAMP MM	12.5 µL	1250 µL
SUBTOTAL VOLUME	23 µL	2300 µL
Sample	2 µL	
REACTION VOLUME	25 µL	

Open EUA's

Provides a path to **establishing generics in diagnostics.**

	 Typical IVD EUA	 CDC EUA	 SalivaDirect™	 SHIELD	 FloodLAMP EasyPCR™	 FloodLAMP QuickColor™
Disclosure of all chemicals and reagents?	No	Yes	Yes	Yes	Yes	Yes
Chemical and reagents available from multiple vendors?	No	Yes	Yes	No	Yes	Yes
Disclosure of primer sequences?	No	Yes std for PCR	Yes	No ProptryThermo	Yes	Yes
Primers commercially available from multiple vendors?	No	Yes	Yes CDC Primers	No ProptryThermo	Yes CDC SD Primers	Yes Available but not launched
Supply chain robust?	No	No/Maybe	Yes	No/Maybe	Yes	Yes
EUA Sponsor Organization Type	For Profit Company	Govt	Academic Not for Profit	Academic Not/For Profit ?	Public Benefit Corp	Public Benefit Corp
Designation of CLIA labs	Kit Sales	N/A open RoR	Impact & Expansion	Impact & Expansion	Impact & Expansion	Impact & Expansion

Push for Regulatory Progress

Prepare for and Respond to Existing Viruses, Emerging New Threats, and Pandemics Act (PREVENT Pandemics Act)

Sec.505.Facilitating the use of real world evidence.	Requires FDA to issue or revise guidance on the use of real-world data and real-world evidence to support regulatory decisionmaking, including with respect to real-world data and real-world evidence from products authorized for emergency use.
Sec.506.Advanced platform technologies.	- Creates an advanced platform technology designation to expedite the development and review of new treatments and countermeasures that use cutting-edge, adaptable platform technologies that can be incorporated or used in more than one drug or biological product. - Requires FDA to issue guidance on the implementation of the new designation.
Sec.507.Increasing EUA decision transparency.	Provides FDA with authority to share more safety and effectiveness information with the public about products authorized for emergency use.
Sec.508.Improving FDA guidance and communication.	- Requires publication of a report identifying best practices across the FDA and other applicable agencies for the development, issuance, and use of guidance documents and for communications with product sponsors and other stakeholders, and a plan for implementing such best practices. - Requires FDA to publish a report on the agency's best practices for communicating with medical product sponsors and other stakeholders, and a plan for implementing such best practices.
Sec.509.GAO study and report on hiring challenges at FDA.	Directs GAO to issue a report assessing FDA's hiring, recruiting, and retention practices, policies, and processes, and their impact on FDA's ability to carry out its public health mission, particularly in light of the COVID-19 pandemic.

Subtitle B—Mitigating Shortages



<https://reaganudall.org/research-funding-opportunity>

Team Effort!



FloodLAMP Team – Thank You!



Randy True
Founder & CEO

- Founder TMI Inc. - Acq. Affymetrix
- Former VP of R&D at Affymetrix

bit.ly/biosketch-randalltrue

co-author global LAMP consortium review: bit.ly/gLAMP-review
<https://abrf.memberclicks.net/jbt-2021-september-issue>



Theresa Ling
UX/Design Lead



Gary Withey, Ph.D.
Vice President of Research & Development



Dr. Peter Antevy, MD
Program Medical Director



Brandon Smith
Lab Assistant

Scientific Advisory Board

Anne Wyllie

Yale, Lead Researcher – SalivaDirect

Bill Hyun

UC San Francisco, Genoa Ventures

Industry Advisors

Tim Lugo

William Blair Biotechnology Group Head

Zarak Khurshid

Asymmetry Capital, Top MDx Analyst

John Edge

Oxford Internet Institute, Blue Field Labs, ID2020

Acknowledgments

Antanas Sadunas, Sam Fogelberg, Carey Tan, Esmé Thornhill-Davis, Jeff Huber, Cliff Wang, Chris Mason, Simon Johnson, Brett Johnson, Mike Finney, Ron Cook

NSVD: Vincent Law, Michael Wells, Katrina Brandis, Gaby Hartley, Tim Earley, Nasreen Haque

Research-Aid Networks: Jeremy Rossman

EMS deployments: Adam, Danny, Alex, Petar, Dr. Paul Pepe

Thanks too to all of our advisors, volunteers & former employees!!

APPENDIX

Successful Commercial Pilots



Commercial deployments in 3 states.



On-site, rapid mass molecular testing performed and operated by EMS staff.



2 cities have opted for mandatory FloodLAMP testing. Prioritizing it over antigen and PCR alternatives.

“

“FloodLAMP’s testing program is a total game-changer. I didn’t know this was possible. I’ve been blown away.”

– *President of Production Company for FloodLAMP Pilot program “ROSA”*


"It's been a godsend!"

– *EMS Training Director running the FloodLAMP program*

”

Collaboration with University of Antioquia (Medellin, Colombia)

Joint NEB /
FloodLAMP project



The New York Times

HEALTH

Several common rapid antigen tests work well for Omicron, according to a new study.

The new findings are from an ongoing U.S. study that began in October and was designed to assess the performance of rapid antigen tests in asymptomatic people.

By **Emily Anthes**

March 1, 2022



medRxiv
THE PREPRINT SERVER FOR HEALTH SCIENCES



BMJ Yale

[View current version of this article](#)

Comparison of Rapid Antigen Tests' Performance between Delta (B.1.61.7;AY.X) and Omicron (B.1.1.529;BA.1) Variants of SARS-CoV-2: Secondary Analysis from a Serial Home Self-Testing Study


 Apurv Soni,  Carly Herbert, Andreas Filippaios, John Broach,  Andres Colubri,  Nisha Fahey, Kelsey Woods, Janvi Nanavati, Colton Wright, Taylor Orwig, Karen Gilliam,  Vik Kheterpal, Thejas Suvarna,  Chris Nowak, Summer Schrader,  Honghuang Lin,  Laurel O'Connor, Caitlin Pretz,  Didem Ayturk,  Elizabeth Orvek,  Julie Flahive, Peter Lazar, Qiming Shi,  Chad Achenbach,  Robert Murphy,  Matthew Robinson, Laura Gibson, Pamela Stamegna,  Nathaniel Hafer,  Katherine Luzuriaga,  Bruce Barton, William Heetderks,  Yukari C. Manabe,  David McManus

doi: <https://doi.org/10.1101/2022.02.27.22271090>

Results From the 7,349 participants enrolled in the parent study, 5,506 met the eligibility criteria for this analysis. A total of 153 participants were RT-PCR+ (61 Delta, 92 Omicron); among this group, 36 (23.5%) tested Ag-RDT+ on the same day and 36 (23.5%) tested Ag-RDT+ within 48 hours as first RT-PCR+. The differences in sensitivity between variants were not statistically significant (same-day: Delta 16.4% [95% CI: 8.2-28.1] vs Omicron 28.2% [95% CI: 19.4-38.6]; and 48-hours: Delta 45.9% [33.1-59.2] vs.

Primer Comparison

All 3 primer sets chosen by HPLC purification



Run with NEB E1700
Fluorimetric LAMP on
QuantStudio7

*Sample is contrived positive with
Gamma BEI spiked into negative
clinical sample*



A = AS1e Rabe-Cepko primer
set for ORF1ab

N = N2 primer set

E = E1 primer set

"H" = HPLC purified

"D" = Desalted