



FloodLAMP QuickColor™ COVID-19 Test Flow




FloodLAMP QuickColor™ COVID-19 Test Diagram

1X Inactivation Saline Solution



1X Inactivation Saline Solution
~50 mL



Amplification Reaction Mix




Reaction Mix Prep: Quick Reference

# Reactions	8	16	24	48
PGS (μL)	92	183	275	550
LAMP MM (μL)	109	218	328	655

+ 23 μL
Reaction Mix
per tube



Pos and Neg
Controls



+ 2 μL
Sample



RESULTS!