


FloodLAMP

Biotechnologies, PBC


Decentralized, Rapid Molecular
Disease Screening for \$1

The Unknowable and How to Prepare for it



David Deutsch, "The Unknowable & how to prepare for it" at TEDxBrussels – ([YouTube](#))

We *still* don't know how to prepare for another variant OR the next pandemic



Ideal Profile of Screening Tests

1

Decentralized

Local control, scalable during surge

2

Rapid

<1hr turn around time

3

MolecularGreater accuracy compared to antigen tests
More quickly adaptable to new targets

4

\$1

Low cost is critical for accessibility and wide distribution

¹ Antigen tests catch only 23% of new asymptomatic infections same day as PCR ; 47% within 48 hrs [Comparison of Rapid Antigen Tests' Performance between Delta \(B.1.61.7; AY.X\) and Omicron \(B.1.1.529; BA1\) Variants of SARS-CoV-2: Secondary Analysis from a Serial Home Self-Testing Study](#)




Programs operating in the real world, at the price of **\$1.25** per person

A Disease Screening Program for Resiliency

1

At-home Accessioned Collection with Mobile App



Streamlines test
processing

2


Family/ Household Pooling




Extends protection

3

Rapid Molecular LAMP Test



<1hr True TAT



Run On or Near Site



LAMP: Loop-Mediated Isothermal Amplification

- Alternative molecular amplification technology
- LAMP generates 10X more DNA product in 1/3 time compared to PCR
- Can be formulated for easy visual/camera-based readout
- Becoming an increasingly important part of diagnostics - many LAMP COVID EUAs!

Recommend [NEB website](https://www.neb.com/) and
global LAMP consortium review:
 (FloodLAMP co-author)

<https://drive.google.com/file/d/1qdf-dL2cmbB7aK4rw57CqjPmWQ9K5X0F/view>

FLOODLAMP ARCHIVE FILE PATH:
 various/glamp/gLAMP Global Consortium - JBT
 Review Paper (Sept 2021).md



	qPCR	LAMP
Equipment	qPCR machine (thermal cycler & optics)	simple heat block or water bath (isothermal)
Time	90 min	30 min
Capital Cost	\$30K-100K	\$300

Molecular Assay Atlas

PEOPLE

SAMPLE

INACTIV

PURIFIC

AMPLIF

READOUT



Individual



Pool: 2-20

Nasopharyngeal
Swab (NP)

Cheek Swab

Anterior Nares
Swab (AN)

Saliva

**TCEP (Chem)
+ Heat**
Rabe Cepko,
HUDSON

**Proteinase K
(Enzyme) + Heat**
Saliva Direct

Heat Only

Commercial
Zymo, Lucigen

Glass Milk
Rabe Cepko

Mag Beads
Kellner, Yu

Skip (Direct)

Commercial
Thermo, Zymo,
Kingfisher, Qiagen

qPCR

LAMP
Colorimetric
Fluorimetric

Misc

RPA, CRISPR, etc.

qPCR Machine
Fluorescence

Visual / Camera
Color Change


Plate Reader
Absorbance
Color Genomics

Lateral Flow Strip




FloodLAMP Biotechnologies,
a Public Benefit Corporation
contact@floodlamp.bio

FloodLAMP Tests



Direct PCR and LAMP Tests



Streamlined Sample Prep

- Shelf stable TCEP/EDTA inactivation solution
- Added directly to dry swabs


- Same sample for both tests
- 2 μ L sample volume



visual readout

QuickColor™ LAMP Test

- Colorimetric LAMP reaction
- 3 primer sets mixed (AS1e, N2, E1)
- 25 minutes at 65°C




Standard RT-qPCR instrument


EasyPCR™ Test

- Multiple PCR Master Mixes Validated
- CDC primers in SalivaDirect config (N1, RNaseP)





Ideal for Austere Sites



FloodLAMP Platform



Scale Up During Omicron Surge



Real World Data



Comparison with diagnostic followup data in progress

- 22 FloodLAMP positives one day in late Dec and only 1 positive on BinaxNOW. All 22 became positive over next few days (reported).
- One site has 100% concordance of FloodLAMP positives with reflex PCR (34/34).

FloodLAMP Surveillance Testing Stats			Thru Mar 2, 2022
Org	Pools	People (Instances)	Positives (incl repeat indiv)
Staff, Preschool, Community	1,107	1,779	44
Conference	61	195	0
Summer Camp	190	696	1
Coral Springs, FL	7,400	21,676	348
Davie, FL	2,235	4,226	55
Bend Fire, OR	617	617	215
TOTAL	11,610	29,189	663



Known adverse results in >10K tests


- 1 suspected false negative (lower instance than discordant lab PCR to FloodLAMP and antigen)
- 2 suspected false positives
- <1% to 5% rerun rate (initial inconclusives that rerun as negative)

Clinical Evaluation Data

Clinical evaluation performed by the Stanford CLIA Lab, with excellent results and praise on the **"really straightforward"** protocol.

EasyPCR™ Test

- 3 copies/μl LoD
- 98% sensitivity (PPA 39/40)
- 100% Specificity (40/40)
- No false positives



Gamma inactivated cell lysate from BEI spiked into raw clinical negative sample

FloodLAMP SwabDirect PCR Result	Comparator Positive	Comparator Negative	Total
Positive	39	0	39
Negative	1	40	41
Invalid	0	0	0
Total	40	40	80
Positive Agreement	97.5% (39/40) 95% CI: 86.8% to 99.9%		
Negative Agreement	100% (40/40) 95% CI: 91.2% to 100%		

QuickColor™ LAMP Test

- 12 copies/μl LoD
- 90% Sensitivity (PPA 36/40)
- Missed positives only high Ct (>36 with direct PCR)
- 100% Specificity (40/40)
- No false positives

FloodLAMP QuickColor Test Result	Comparator Positive	Comparator Negative	Total
Positive	36	0	36
Negative	4	40	44
Total	40	40	80
Positive Agreement	90.0% (36/40) 95% CI: 76.3% to 97.2%		
Negative Agreement	100% (40/40) 95% CI: 91.2% to 100%		

Source of Specimens: Stanford COVID-19 Clinical Testing Program
 Specimen Type: Anterior Nares Swab in PBS, previously tested and frozen
 Comparator Test: Hologic Panther Fusion SARS-CoV-2 Assay and Hologic Panther Aptima SARS-CoV-2 Assay

EUA Submissions – 2 Full + Pre EUA

FloodLAMP QuickColor™ COVID-19 Test

Instructions for Use v1.2

IVD

COVID-19 Emergency Use Authorization Only
For *in vitro* diagnostic (IVD) Use

www.floodlamp.bio
FloodLAMP Biotechnologies, PBC | 930 Brittan Ave. San Carlos, CA 94070 USA

FloodLAMP EasyPCR™ COVID-19 Test

Instructions for Use v1.1

IVD

COVID-19 Emergency Use Authorization Only
For *in vitro* diagnostic (IVD) Use

www.floodlamp.bio
FloodLAMP Biotechnologies, PBC | 930 Brittan Ave. San Carlos, CA 94070 USA



**FloodLAMP
Biotechnologies**

A Public Benefit Corporation

**FloodLAMP Pooled Swab
Collection Kit DTC**

For use with the FloodLAMP
Mobile App

Open EUA – Disclosure of Components

In **contrast** to industry standard trade secret primer sequences and reaction components.

Table 1: Validated reagents used with the Test

Item	Concentration	Chemical Composition	Vendor	Catalog Number
TCEP	.5 M	tris(2-carboxyethyl)phosphine hydrochloride	Sigma-Aldrich	646547-10XIML
EDTA	.5 M	Ethylenediaminetetraacetic acid	Thermo Fisher	15575020
NaOH	10 N	Sodium Hydroxide	Sigma-Aldrich	SX0607N-6
Nuclease-free Water		Ultrapure Water, nuclease-free	Thermo Fisher	10977015
NaCl	5 M	Sodium Chloride	Thermo Fisher	24740011
Guanidine HCl	6 M	Guanidine Hydrochloride	Sigma-Aldrich	SRE0066
Colorimetric LAMP MM*		Colorimetric LAMP Master Mix	New England Biolabs	M1804

* Item may not be substituted for equivalent. Only the specified vendor and catalog number may be utilized.

Table 2: Primer names and sequences

Primer Name	Sequence (5'–3')
ORFlab gene (AS1e)	
Orflab_FIP	TCAGCACACAAAGCCAAAAATTTATTTTCTGTGCAAAAGGAAATTAAGGAG
Orflab_BIP	TATTGGTGGAGCTAAACTTAAAGCCTTTTCTGTACAATCCCTTTGAGTG
Orflab_F3	CGGTGGACAAATTGTCAC
Orflab_B3	CTTCTCTGGATTAAACACACTT
Orflab_LF	TTACAAGCTTAAAGAATGTCTGAACACT
Orflab_LB	TTGAATTTAGGTGAAACATTGTCACG
N Gene (N2)	
N2_FIP	TTCCGAAGAACGCTGAAGCGGAACCTGATTACAAACATTGGCC
N2_BIP	CGCATTGGCATGGAAGTCACAATTTGATGGCACCTGTGTA
N2_F3	ACCAGGAACATAACAGACAAG
N2_B3	GACCTGATCTTTGAAATTTGGATCT
N2_LF	GGGGGCAAAATGTGCAATTTG
N2_LB	CTTCGGGAACGTGGTTGACC

Table 7: Primer-Guanidine Solution

Component	Volume (1 reaction)	Volume (1 reaction x 100) 1 x 96-plate w/ 4% overage
10X LAMP Primer Mix	2.5 µL	250 µL
Guanidine HCl (400 mM)	2.5 µL	
Guanidine HCl (6 M)		16.7 µL
Nuclease-free Water	5.5 µL	783 µL
TOTAL VOLUME	10.5 µL	1050 µL

Table 8: Colorimetric LAMP Amplification Reaction


Component	Volume (1 reaction)	Volume (100 reactions)
Primer-Guanidine Solution	10.5 µL	1050 µL
Colorimetric LAMP MM	12.5 µL	1250 µL
SUBTOTAL VOLUME	23 µL	2300 µL
Sample	2 µL	
REACTION VOLUME	25 µL	

Open EUAs

Provides a path to **establishing generics in diagnostics**.

Enables **global dissemination** of highly respected FDA authorized tests.

	 Typical IVD EUA	 CDC EUA	 SalivaDirect™	 SHIELD	 FloodLAMP EasyPCR™	 FloodLAMP QuickColor™
Disclosure of all chemicals and reagents?	No	Yes	Yes	Yes	Yes	Yes
Chemical and reagents available from multiple vendors?	No	Yes	Yes	No	Yes	Yes
Disclosure of primer sequences?	No	Yes std for PCR	Yes	No ProptryThermo	Yes	Yes
Primers commercially available from multiple vendors?	No	Yes	Yes CDC Primers	No ProptryThermo	Yes CDC SD Primers	Yes Available but not launched
Supply chain robust?	No	No/Maybe	Yes	No/Maybe	Yes	Yes
EUA Sponsor Organization Type	For Profit Company	Govt	Academic Not for Profit	Academic Not/For Profit ?	Public Benefit Corp	Public Benefit Corp
Designation of CLIA labs	Kit Sales	N/A open RoR	Impact & Expansion	Impact & Expansion	Impact & Expansion	Impact & Expansion



FloodLAMP

Biotechnologies, PBC

Decentralized, Rapid, Molecular,
Disease Screening for \$1